

FULL IMPLANTATION SURGICAL KIT BOX

This Full Implantation Kit is intended for convenient access to surgical instruments during bone preparation and implant placement procedures, inside the Kit you'll find all the instruments required for a wide range of implant placement cases.

The Kit is comfortable to use, intuitively organized, and clearly labeled, making it easy to put them back in the tray after use.

Please follow the recommended drilling protocol to ensure optimal outcome.

FEATURES:

- Drills are Carbon DLC coated, enabling reduced resistance and heating during drilling, for a more secure, precise, and stable drill.
- Removable Tray, making it easy to clean and maintain.
- Suitable for autoclaving, ensuring proper sterilization - at a temperature of 134°C.
- All instruments are clearly labeled, making them easy to identify and use.
- The tools are made from durable Stainless Steel, ensuring high-quality performance.



KIT CONTENT			CODE: SK006	
ITEM	NAME	L (mm)	D (mm)	CODE
	Marking Drill	16	Ø1.9	MDL19
	Lance Drill	16	Ø1.5	LPD15
	Drill Extender	29	-----	TDE
	Countersink Drill	16	Ø3.75-Ø4.2	CD3742
	Lindemann Drill	29	Ø5.0-Ø6.0	LD1520
	Round Ball Bur	30	Ø3.8	RB38
	Parallel Pin	10	-----	PP10
	Parallel Pin	16	-----	PP16
	Standard Drill	16	Ø2.0	TD200
			Ø2.8	TD280
			Ø3.2	TD320
			Ø3.65	TD365
			Ø4.0	TD400
			Ø4.2	TD420
			Ø5.2	TD520
	Hex Driver	7	Ø1.25	HD1207
	Hex Driver	15	Ø1.25	HD1215
	Hex Driver	15	Ø2.0	HD2015-SL
	Hex Driver	10	Ø2.42	HD2410
	Hex Driver	15	Ø2.42	HD2415
	Motor Mount	20	Ø1.25	MM1220
	Motor Mount	28	Ø1.25	MM1228
	Motor Mount	28	Ø2.0	MM2028-SL
	Motor Mount	20	Ø2.42	MM2420
	Motor Mount	28	Ø2.42	MM2428
	Hand Hex Driver	7	Ø1.25	HHD1207
			Ø1.25	HHD1215
	Integral Step&Stop Drill		Ø2.0-Ø2.8	SD2028L6SC
			Ø2.0-Ø2.8	SD2028L8SC
			Ø2.0-Ø2.8	SD2028L10SC
			Ø2.0-Ø2.8	SD2028L11SC
			Ø2.0-Ø2.8	SD2028L13SC
			Ø2.8-Ø3.2	SD2832L6SC
			Ø2.8-Ø3.2	SD2832L8SC
			Ø2.8-Ø3.2	SD2832L10SC
			Ø2.8-Ø3.2	SD2832L11SC
			Ø2.8-Ø3.2	SD2832L13SC
			Ø3.2-Ø3.65	SD32365L6SC
			Ø3.2-Ø3.65	SD32365L8SC
			Ø3.2-Ø3.65	SD32365L10SC
	Ø3.2-Ø3.65	SD32365L11SC		
	Ø3.2-Ø3.65	SD32365L13SC		
	Ø3.65-Ø4.0	SD36540L6SC		
	Ø3.65-Ø4.0	SD36540L8SC		
	Ø3.65-Ø4.0	SD36540L10SC		
	Ø3.65-Ø4.0	SD36540L11SC		
	Ø3.65-Ø4.0	SD36540L13SC		
	Ø4.0-Ø4.5	SD4045L6SC		
	Ø4.0-Ø4.5	SD4045L8SC		
	Ø4.0-Ø4.5	SD4045L10SC		
	Ø4.0-Ø4.5	SD4045L11SC		
	Ø4.0-Ø4.5	SD4045L13SC		
	Torque Wrench	-----	Ø6.35	STRA
				TRA
	Depth Probe	-----	-----	DEGG

25-100 Ncm
10-45 Ncm

• Product image size and proportions are for illustration purposes only. Please refer to the table for actual sizes and proportions.